



as part of this project explain what needs to be done on day 9 - provide easy breakdown of what needs to be done as I am a shadow Product manager

Day 9 focuses on shipping a minimal end-to-end UI in Streamlit with simple authentication so the demo can run from note input to structured summary with basic role-based views.^[1]

Core outcome

A working Streamlit app with three tabs (Upload/Enter Note, Summarize, Logs) and simple, file-based auth using streamlit-authenticator with clinician/admin roles, wired to the existing de-id → index → retrieve → summarize pipeline artifacts from Days 5–8.^[1]

What to build

- Streamlit app file structure with app/main.py containing 3 tabs for input, summarization, and viewing logs to demonstrate the full workflow as a product demo entry point.^[1]
- Authentication via streamlit-authenticator using a small YAML config with hashed passwords and role labels for clinician and admin to enable basic access restrictions suitable for a portfolio MVP.^[1]
- Wiring to existing modules: deidentification, indexer, retriever, summarizer so the UI triggers these in sequence or skips already-indexed notes, then offers a downloadable summary artifact.^[1]

Step-by-step tasks

- Create Streamlit scaffolding: set page title, sidebar login, and three tabs named Upload/Enter Note, Summarize, and Logs to mirror the defined Day 9 objectives and artifacts.^[1]
- Add streamlit-authenticator setup: prepare a config.yaml with usernames, bcrypt-hashed passwords, and roles (clinician, admin), then instantiate the authenticator and gate the app's content behind a successful login per the plan.^[1]
- Implement Upload/Enter Note tab: allow pasting text or uploading a .txt file; provide a button to run de-identification and either index immediately or let the user choose to skip if already indexed, consistent with the Day 9 instruction to deid/index or skip if already indexed.^[1]

- Implement Summarize tab: accept a selected note or free text, invoke the RAG retrieval from Day 7 and the prompt-based summarizer from Day 8, and show a structured clinical summary, then enable download as a file for portfolio demos as indicated for downstream polish steps.^[1]
- Implement Logs tab: display recent app actions for demo transparency; for Day 9 this can be a simple in-app log view placeholder prior to the hash-chained audit in Day 10, aligning with the “Logs” tab requirement.^[1]

Role checks

- Clinician role: can upload/enter notes, trigger de-identification, indexing, and generate summaries within the UI per the primary user path described in the earlier weeks.^[1]
- Admin role: can additionally view the Logs tab content to prepare for Day 10’s audit logging and administrative views, as hinted by the role split in Day 9 and expanded in Days 11–12.^[1]

Files to produce

- app/main.py implementing Streamlit tabs, login, and calls into the existing pipelines to meet the Day 9 demo goal.^[1]
- app/streamlit_config.yaml containing users, hashed passwords, and role mapping used by streamlit-authenticator as specified in Day 9.^[1]

Demo flow to validate

- Log in as clinician using the YAML-configured account, paste a synthetic note, click to de-identify and index, then navigate to Summarize and generate a structured summary that adheres to the Day 8 prompt format for clinical sections.^[1]
- Log in as admin and confirm access to the Logs tab view to verify basic role separation in preparation for later security and audit tasks scheduled in Week 3 and Day 10.^[1]

✱

tell me in easy step by step breakdown how to achieve above requirements, i dont know how to do it

Great—let’s make Day 9 doable. Your goal: a simple Streamlit app with login and three tabs wired to your existing de-id → index → retrieve → summarize pipeline. We’ll move in small steps and you can paste back errors as they appear.

Step 0 — Pre-checks (5–10 min)

- Confirm your repo folders: docs/ and app/ exist.
- Ensure your Day 5–8 files are present/importable: app/deid_pipeline.py, app/indexer.py, app/rag_pipeline.py, app/summarizer.py.

- Create/activate your Python venv. Then: pip install streamlit streamlit-authenticator pyyaml bcrypt.

Step 1 — Minimal Streamlit scaffold (15 min)

- In app/, create main.py with a blank Streamlit page that runs:

```
# app/main.py
import streamlit as st
st.set_page_config(page_title="Clinical Summarizer", layout="wide")
st.title("HIPAA-aware RAG Summarizer (MVP)")
```

- Run: streamlit run app/main.py
- You should see a blank page with the title. If this works, proceed.

Step 2 — Add simple auth with roles (20–30 min)

- In app/, create streamlit_config.yaml with two users and roles:

```
credentials:
  usernames:
    clinician1:
      email: clinician1@example.com
      name: Clinician One
      password: "$2b$12$replace_with_bcrypt_hash"
      role: clinician
    admin1:
      email: admin1@example.com
      name: Admin One
      password: "$2b$12$replace_with_bcrypt_hash"
      role: admin
  cookie:
    expiry_days: 1
    key: some_random_secret
    name: auth_cookie
  preauthorized:
    emails: []
```

- Generate bcrypt hashes (don't store plaintext). In a Python shell:

```
import bcrypt
print(bcrypt.hashpw(b"your_password", bcrypt.gensalt()).decode())
```

- Wire authenticator in main.py:

```
import yaml
from yaml.loader import SafeLoader
import streamlit_authenticator as stauth

with open("app/streamlit_config.yaml") as f:
    config = yaml.load(f, Loader=SafeLoader)

authenticator = stauth.Authenticate(
```

```

        config["credentials"],
        config["cookie"]["name"],
        config["cookie"]["key"],
        config["cookie"]["expiry_days"],
    )

    name, auth_status, username = authenticator.login(location="sidebar")

    if auth_status is False:
        st.error("Invalid username or password")
    elif auth_status is None:
        st.info("Please log in")
    else:
        st.sidebar.success(f"Logged in as {name}")
        role = config["credentials"]["usernames"][username]["role"]
        st.session_state["role"] = role
        authenticator.logout("Logout", location="sidebar")
        # show tabs here (Step 3)

```

- Test login with your hashed creds. If you can log in and see the sidebar status, continue.

Step 3 — Create three tabs (10–15 min)

- After successful login block, add tabs:

```
upload_tab, summarize_tab, logs_tab = st.tabs(["Upload/Enter Note", "Summarize", "Logs"])
```

- We'll fill each tab in Steps 4–6.
- Role gating: only admins see Logs.

```
show_logs = st.session_state.get("role") == "admin"
```

Step 4 — Upload/Enter Note tab: paste or upload, then de-identify + index (30–45 min)

- Inside upload_tab:

```

with upload_tab:
    st.subheader("Enter or Upload Note")
    note_text = st.text_area("Paste note text", height=200)
    file = st.file_uploader("...or upload a .txt file", type=["txt"])
    if file and not note_text:
        note_text = file.read().decode("utf-8", errors="ignore")

    col1, col2 = st.columns(2)
    with col1:
        run_deid = st.button("De-identify & Index")
    with col2:
        skip_index = st.button("Skip (already indexed)")

```

- Wire to your de-id and indexer modules. Adjust import names to match your files.

```

from deid_pipeline import deidentify_note # returns deid_text, encrypted_span_map_path
from indexer import index_note # takes deid_text, returns note_id

if run_deid and note_text:
    with st.spinner("De-identifying..."):
        deid_text, enc_map_path = deidentify_note(note_text)
        st.success("De-identified.")
        st.text_area("De-identified preview", deid_text, height=160)

    with st.spinner("Indexing..."):
        note_id = index_note(deid_text)
        st.session_state["last_note_id"] = note_id
        st.session_state["last_deid_text"] = deid_text
        st.success(f"Indexed note_id: {note_id}")

elif skip_index and note_text:
    # If user knows it was indexed earlier, just store text for summarize step
    st.session_state["last_deid_text"] = note_text
    st.info("Skipped indexing; text saved for summarization.")

```

- If your functions have different names/signatures, adapt accordingly.

Step 5 — Summarize tab: retrieve + summarize, then download (30–45 min)

- Inside summarize_tab:

```

from rag_pipeline import retrieve_context # input: text or note_id → list[str]
from summarizer import generate_summary # input: context → structured summary str

with summarize_tab:
    st.subheader("Summarize")
    source_choice = st.radio("Use:", ["Last de-identified text", "Note ID"], horizontal=True)
    user_note_id = st.text_input("Note ID", value=st.session_state.get("last_note_id", ""))

    if st.button("Generate Summary"):
        with st.spinner("Retrieving context..."):
            if source_choice == "Note ID" and user_note_id:
                contexts = retrieve_context(note_id=user_note_id)
            else:
                deid_text = st.session_state.get("last_deid_text", "")
                contexts = retrieve_context(text=deid_text)
        with st.spinner("Summarizing..."):
            summary = generate_summary(contexts)
            st.success("Summary ready")
            st.text_area("Structured Summary", summary, height=250)
            st.download_button("Download .txt", data=summary, file_name="summary.txt")

```

- Keep output strictly from context per your Day 8 prompt design.

Step 6 — Logs tab: simple placeholder, admin-only (10–15 min)

- In logs_tab, show recent events from a simple JSONL file if available; otherwise a placeholder. Only show if role == admin:

```

with logs_tab:
    if st.session_state.get("role") != "admin":
        st.info("Admins only.")
    else:
        st.subheader("Logs (preview)")
        try:
            with open("logs/app_audit.jsonl") as f:
                lines = f.readlines()[-100:]
                for line in lines:
                    st.code(line)
        except FileNotFoundError:
            st.write("No logs yet. Day 10 will add hash-chained audit logs.")

```

- Day 10 will replace this with hash-chained logs.

Step 7 — Minimal UX polish and state glue (10–15 min)

- After actions, show clear success/info messages.
- Store `last_note_id` and `last_deid_text` in `st.session_state` so users can move from Upload → Summarize smoothly.
- Add a small help expander on each tab explaining the step.

Step 8 — Smoke test the demo flow (10–20 min)

- Login as clinician → paste a synthetic note → De-identify & Index → Summarize → download.
- Login as admin → confirm Logs tab appears (even if only a placeholder today).

Common pitfalls and quick fixes

- Module import errors: ensure `main.py` and your pipeline files are importable (`PYTHONPATH` includes `app/`). Quick fix: run from repo root and use `python -m streamlit run app/main.py`.
- `bcrypt` hash not accepted: confirm you pasted the full hash (starts with `2b`). Regenerate if unsure.
- File uploads reading empty: ensure `file.read()` and `decode` call are used once; reset the uploader if needed.
- State lost on rerun: write to `st.session_state` right after each action and check the keys before using them.

Your turn

- Which step do you want to attempt first—auth setup or the Upload tab wiring? Tell me where you are, paste any error you hit, and we'll fix it together one line at a time.

✱✱

1. COMPLETE-WEEK-BY-WEEK-TASKS.docx
2. COMPLETE-WEEK-BY-WEEK-TASKS.docx

